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

Hand-Held Refractometer Operating Instructions

MODEL : BR32T、BR62T、BR90T、SA10T、SA28T、AL80、SO25、CF18

I. Part

1. Prism
2. Cover plate
3. Zero correct screw
4. Mirror tube
5. Focus adjuster

II. OPERATING PROCEDURE

1. Take the prism (1) end of refractometer towards the bright light and turn focus adjuster (5) clockwise or counterclockwise until the reticule can be seen clearly.

2. Zero adjustment

- for **BR32、SA10T、SA28T、AL80、SO25、CF18**

Lift the cover plate (2) up, drop 1 ~ 2 drops of pure water on the prism. Cover the plate and press it lightly and make sure there's no bubbles. Adjusts the zero correct screw (3) to match zero boundaries between bright and dark side coincide with the null line.

- for **BR62T**

Replace pure water with saturated sodium chloride solution as zero standard. The refractive indexes of the saturated sodium chloride solution which refractive indexes are 29.9% (15°C), 29.6% (20°C), 29.2% (25°C).

- for **BR90T**

Drop 1 drop of dioptric oil on the bright and matte surface of the reference block. Lift the cover plate (2) up, put the bright side of reference block on the surface of the prism and fix it. (Do NOT cover the cover plate.) Turn and adjust the zero correct screw (3) to have the bright/dark boundary located at the Brix 78.8% reference line.

3. Measurement: Lift the cover plate (2) up, clean surface of prism with tissue paper, then drop 1~2 drops of solution for measurement and make sure there's no bubbles. Cover the plate, press it lightly, then read the corresponding scale, where the bright/dark boundary located, scale reading brix of sample.
4. After measurement, clean and wipe the surface of prism and cover plate with wet tissue paper. Please store the refractometer well after it is dry.

III. MAINTENANCE

1. The temperature of zeroing standard should be the same with ambient temperature. If ambient temperature varies greatly, zeroing should be done again.
2. Never flush the refractometer, to prevent water from draining into the pipe of refractometer.
3. It should be kept in the dry, clean and non-corrosiveness environment.
4. Do not touch and scratch optical surfaces.
5. Do not shock strongly during transportation.

※ **TEMPERATURE COMPENSATION:** The measuring results will be affected slightly when temperature is not at 20°C, temperature correction table listed below used for temperature compensation.

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