

SS, Stainless Steel Filter Holder



Features


- Inner spin-lock connection (M381450) facilitates fast and stable installation without clamp.
- Intuitive operation by 360° ring graduation mark.
- SS316 funnel offers excellent chemical resistance and can be fast sterilized by flame.

Specifications

Material	100 ml	300 ml	500 ml
Funnel	SS316	SS316	SS316
Support base	SS316	SS316	SS316
Membrane support	SS316	SS316	SS316
Funnel lid	SS316	SS316	SS316

Specifications

Funnel capacity	100 ml	300 ml	500 ml
Filter diameter	47 mm	47 mm	47 mm
Effective filtration area	9.6 cm ²	9.6 cm ²	9.6 cm ²
Support pore size	70 μm	70 μm	70 μm



Ordering Information

SS 1, Stainless Steel Filter Holder (100 ml) 167103-10

1	Stainless steel funnel, 100 ml	167103-61
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.8)	167110-16

SS 3, Stainless Steel Filter Holder (300 ml) 167103-30

1	Stainless steel funnel, 300 ml	167103-63
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.8)	167110-16

SS 5, Stainless Steel Filter Holder (500 ml) 167103-50

1	Stainless steel funnel, 500 ml	167103-65
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.8)	167110-16

SS, Stainless Steel Lid

100 ml	Stainless steel lid	167103-11
300 ml	Stainless steel lid	167103-13
500 ml	Stainless steel lid	167103-15