


操作說明書

Cat. No. 197000-60

A. 簡介

抽氣血清瓶蓋組為一相容於 GL45 規格血清瓶之蓋體，可將實驗室常用之血清瓶轉接作為廢液接收之用。組裝簡單且可高溫滅菌重複使用，但請勿以明火直接燒烤。

B. 結構說明



** 參照Fig.2，液體端(LIQ)接抽吸器，氣體端(AIR)接真空源即可使用。

C. 訂購資訊

197000-60	抽氣血清瓶蓋組
197000-60-500	500ml抽氣血清瓶組
197000-63	血清瓶蓋矽膠墊環
199100-68-5	矽膠管 (Φ4.2xΦ8mm), 5公分
197100-01	Lafil 100 生化廢液抽吸器
197405-11(22)	Lafil 400 Plus 生化廢液抽吸器
199100-00	BioDolphin 廢液抽吸套件組(淺藍)


INSTRUCTION MANUAL

Cat. No. 197000-60

A. Introduction

The GL45 suction cap storage bottle can make a typical tissue bottle work as a suction bottle. It is autoclavable, reusable, easily assembled and ideal for general waste liquid suction.

B. Diagram



** See Fig.2: Connect “LIQ” (Liquid inlet) to suction kit and connect “AIR” (air outlet) to vacuum source.

C. Order Information

197000-60	GL45 suction cap storage bottle
197000-60-500	500ml GL45 storage bottle with suction cap
197000-63	Silicone O-ring for GL45 storage bottle
199100-68-5	Silicone tube (Φ4.2xΦ8mm), 5cm
197100-01	Lafil 100 portable suction system
197405-11(22)	Lafil 400 Plus portable suction system
199100-00	BioDolphin suction kit (light blue color)