

Large Volume Suction System

All in One · Safe · Space-Saving



Features



Safe Operation & Overflow Protection

- 1 ♦ The fence-like platform prevents negligence from tipping over the waste bottle.
- ♦ Prote™ overflow protection prevents liquid from being sucked into the pump.



All in One & Space-Saving

- 2 ♦ Integrates vacuum source with waste bottle and suction kit, saves 50% space.
- ♦ Allows 2 suction kits to operate at the same time.



Large Waste Bottle

- 3 ♦ 4L PP heavy duty vacuum bottle with quick-coupling connectors
- ♦ 3.5L laboratory glass bottle (optional)



Certification

CE / RoHS / IP 41



▲ 3.5L laboratory glass bottle (Optional)

Applications

- ♦ Cell culture aspiration
- ♦ Laboratory waste suction

Specification

Model	Lafil 200 - Plus / BioDolphin	Lafil 200 ^{eco} - Plus / BioDolphin
Power	DC12V, 1.5A	
Vacuum range	1013 ~ 213 mbar	880 ~ 370 mbar
Max. flow rate	10 L/min	
Aspiration rate	17 mL/sec	
Waste bottle	• 4L PP vacuum bottle with quick-coupling connectors • 3.5L laboratory glass bottle (optional)	
Auto shut-off	-	Yes
Noise level	55 dB	
Weight	2.1 Kg	2.3 Kg
Dimensions (LxWxH)	37.5 x 25.8 x 21.7 cm	
	Plus	BioDolphin
Suction kit	• Handle • 1-channel / 8-channel tip adaptor with ejector • 1-channel tip adaptor (150 mm)	• Handle • 1-channel / 8-channel tip adaptor with ejector • 1-channel tip adaptor (50 / 150 mm) • 1-channel / 8-channel SS adaptor (40 mm) • 1-channel SS adaptor (80 mm) • Rack



Lafil 200 - Plus, Suction System



Lafil 200 - BioDolphin, Suction System



**Ordering
Information**

197205-01(02)	Lafil 200-Plus, Suction System with AC100~240V adaptor, US plug (EU plug)
197203-01(02)	Lafil 200-BioDolphin, Suction System with AC100~240V adaptor, US plug (EU plug)
197215-01(02)	Lafil 200 ^{eco} -Plus, Suction System with AC100~240V adaptor, US plug (EU plug)
197213-01(02)	Lafil 200 ^{eco} -BioDolphins, Suction System with AC100~240V adaptor, US plug (EU plug)
197200-53	4L PP Vacuum Bottle with Quick-Coupling Connectors
197000-60-3500	3.5L Laboratory Glass Bottle with GL45 Suction Cap
197000-60	Suction Cap for GL45 Bottle
167200-41	Universal Foot Switch (Momentary type)
167200-41-1	Universal Foot Switch (Maintained type)

Rocker Scientific Co., Ltd.

Headquarter

19F.-6, No.206, Guanghua 1st Rd., Lingya Dist.,
Kaohsiung City 802, Taiwan (R.O.C.)
T/ +886-7-227-3358 F/ +886-7-227-3716

Taipei office

11F., No.6, Zhulin Rd., Linkou Dist.,
New Taipei City 244, Taiwan (R.O.C.)
T/ +886-2-2603-3311 F/ +886-2-2603-6622

Get more information for Rocker

www.rocker.com.tw
export@rocker.com.tw

