



WaterVac 200
Vacuum Filtration System

Efficient / Space-Saving / Versatile

Compliance with ISO 8199

Features



Save Time

- ◆ Integrates liquid pump and filter holders into an instrument, transferring filtered liquid directly to waste.
- ◆ Proceeds 2 tests simultaneously.



Save 70% Bench Space

- ◆ No need of waste bottle and tubings saving much bench space.



Versatile

- ◆ Replaceable funnel adaptor makes it compatible with various kinds of filter holders.



Certification

CE / RoHS / IP 41



▲ Compatible with various kinds of filter holders.

Applications

- ◆ Microbiological test by Membrane Filter Technique
- ◆ Suspended solid test for water monitoring

Specification

Model	WaterVac 200		
Power	DC24V, 1.8A		
Max. vacuum	313 mbar		
Max. liquid flow rate	4.5 L/min		
Motor	60W Brushless DC Motor, improves the pump life and performance		
Noise level	60 dB		
Hose barb	ID8		
Material contact with filtrate	Stainless steel, FKM, Silicone, PP, PPS		
International certification	CE, RoHS, IP 41		
Dimension (LxWxH)	21.4 x 18.3 x 26.4 cm		
	WaterVac 200 - MS	WaterVac 200 - MB	WaterVac 201 - MB
Adaptor	Stopper	Spin-Lock	Spin-Lock
Filter holder / funnel	MF 3 / MF 5, Magnetic Filter Holder (Optional)	SF Stainless Steel Funnel (Optional)	SF Stainless Steel Funnel (100 ml x 2)



WaterVac 200 - MS
Vacuum Filtration System



WaterVac 200 - MB
Vacuum Filtration System



WaterVac 201 - MB
Vacuum Filtration System



Ordering Information

190210-01(02)	WaterVac 200-MS, Vacuum Filtration System with No. 8 Stopper Funnel Adaptor, AC100~240V Adaptor, US Plug (EU Plug)
190220-01(02)	WaterVac 200-MB, Vacuum Filtration System with Rocker SF Funnel Adaptor, AC100~240V Adaptor, US Plug (EU Plug)
190221-01(02)	WaterVac 201-MB, Vacuum Filtration System with Rocker SF Funnel Adaptor, SF Stainless Steel Funnel 100 ml x 2 pcs, AC100~240V Adaptor, US Plug (EU Plug)
180100-01	SF Stainless Steel Funnel 100 ml
180100-03	SF Stainless Steel Funnel 300 ml
180100-05	SF Stainless Steel Funnel 500 ml
200300-01	MF 3, Magnetic Filter Holder 300 ml
200500-01	MF 5, Magnetic Filter Holder 500 ml
180300-30	SS316 Adaptor for Pall disposable funnel
190100-73	SS316 Adaptor for Millipore disposable plastic funnel
190100-74	SS316 Adaptor for Sartorius disposable plastic funnel

Rocker Scientific Co., Ltd.

Headquarter

19F.-6, No.206, Guanghua 1st Rd., Lingya Dist.,
Kaohsiung City 802, Taiwan (R.O.C.)
T/ +886-7-227-3358 F/ +886-7-227-3716

Taipei office

11F., No.6, Zhulin Rd., Linkou Dist.,
New Taipei City 244, Taiwan (R.O.C.)
T/ +886-2-2603-3311 F/ +886-2-2603-6622

Get more information for Rocker

www.rocker.com.tw
export@rocker.com.tw

