

# STREAM

## Continuous Filtration System

Effort-Saving • Quick • User-Friendly



### Features



#### Continuous Filtration

- 1
  - ◆ Ensure an unceasing process of large capacity solution.
  - ◆ Anchor design allows the maximum suction to the last bit of liquid.



#### Effort-Saving & User-Friendly

- 2
  - ◆ No need of transferring or pouring solution manually.
  - ◆ Completes large volume filtration in one set, no more redundant work.



#### Chemical Resistance

- 3
  - ◆ Glass and PTFE construction ensures highly chemical resistance.
  - ◆ Autoclavable components avoid the possibility of contamination.



#### STREAM

Continuous Filter Holder includes:

- (1) Continuous Filtration Kit
- (2) 300 ml or 1000 ml Filter Holder
- (3) GL45 Filtration Adaptor

### Applications

- ◆ Solvent purification
- ◆ Large volume filtration


## Specification

Model	Stream 3	Stream 10
Funnel capacity	300 ml	1000 ml
Filter diameter	47 mm	90 mm
Hose barb	ID4~ID8	ID4~ID8
Part	Materia	Material
Continuous filtration kit	PTFE, Silicone	PTFE, Silicone
Filter holder	Borosilicate glass	Borosilicate glass
Membrane support	SS316	SS316
Clamp	Aluminum	Aluminum
Filtration adaptor	PP	PP

## Examples

Compatible with

◆ Various GL45 Laboratory Glass Bottle



Stream 3 / 10

### Stream 3, Continuous Filter Holder (300 ml) 167200-S3

1	Continuous Filtration Kit for 300 ml Glass Filter Holder	167210-36
2	300 ml Glass Filter Holder	167120-33
3	Filtration Adaptor for GL45 Bottle	197000-65
optional	Sintered Glass Support Base (Including No.8 Silicone Stopper)	167230-03

### Stream 10, Continuous Filter Holder (1000 ml) 167200-S10

1	Continuous Filtration Kit for 1000 ml Glass Filter Holder	167210-35
2	1000 ml Glass Filter Holder	167120-37
3	Filtration Adaptor for GL45 Bottle	197000-65
optional	Sintered Glass Support Base (Including No.8 Silicone Stopper)	167230-11

<b>167300-11(22)</b>	Rocker 300, Vacuum Pump, AC110V, 60Hz (AC220V, 50Hz)
<b>167400-11(22)</b>	Rocker 400, Vacuum Pump, AC110V, 60Hz (AC220V, 50Hz)
<b>189300-11(22)</b>	Rocker 300C, PTFE Coated Chemical Resistant Vacuum Pump, AC110V, 60Hz (AC220V, 50Hz)
<b>189400-11(22)</b>	Rocker 400C, PTFE Coated Chemical Resistant Vacuum Pump, AC110V, 60Hz (AC220V, 50Hz)
<b>197000-55</b>	Laboratory Glass Bottle, 5000 ml



## Rocker Scientific Co., Ltd.

### Headquarter

19F.-6, No.206, Guanghua 1st Rd., Lingya Dist.,  
Kaohsiung City 802, Taiwan (R.O.C.)  
T/ +886-7-227-3358 F/ +886-7-227-3716

### Taipei office

11F., No. 402, Sec. 1, Ren'ai Rd., Linkou Dist.,  
New Taipei City 244014, Taiwan (R.O.C.)  
T/ +886-2-2603-3311 F/ +886-2-2603-6622

Get more information for Rocker

[www.rocker.com.tw](http://www.rocker.com.tw)  
[export@rocker.com.tw](mailto:export@rocker.com.tw)

