

Oil Free Vacuum Pump / Compressor



Features

► 2-in-1 machine

Rocker 430 works as both vacuum pump and compressor. Equipped with gauges and regulators for easy monitoring and adjusting the vacuum level or delivery pressure.

► Oil-free & care-free

Rocker pumps are driven by piston that needs no lubricant, making them maintenance-free and oil-mist-free.

► Prote™ overflow protection

Equipped with Prote™ overflow protection to prevent liquid from being sucked into the pump.

► Overheat protection

All Rocker pumps are equipped with overheat protection device that shuts down when overheated and resumes once cooled down, making the pumps safer and long-lasting.

► Moisture trap with filter cartridge

Filter cartridge filters particles and moisture to ensure air cleanness of both inlet and supply.

► International certification

CE certification
IP 30 protection
RoHS certification

► Two-Year Warranty

Applications

- Vacuum filtration
- Microbiological test
- Suspended solid test
- Air sampling
- Pressure filtration

Ordering Information

167430-11(22)

Rocker 430, Vacuum Pump
AC110V, 60Hz (AC220V, 50Hz)



Rocker 430

Specification

Model	Rocker 430	
PERFORMANCE DATA		
Max. vacuum	145 mbar abs.	
Max. pressure	80 psig	
Max. flow rate	34 L/min @60Hz	28 L/min @50Hz
Hose barb	ID8	
Noise level	52 dB	
Net weight	5.5 Kg	
Dimension (LxWxH)	25 x 24 x 20 cm	
ELECTRICAL DATA		
Voltage	110V	220V
Frequency	50 / 60Hz	
Max. power	147W (1/6 HP)	125W (1/6 HP)
Max. current	1.5A	0.6A
Safety	Overheat Protection	
OTHER DATA		
Equipped with	Vacuum regulator / Pressure regulator Overload Protection / Filter cartridge	