

# Oil Free Compressor



WATCH VIDEO

2  
YEAR

WARRANTY

## Features

### ► Compact & care-free

Rocker oil-free compressors are compact in size, light in weight and oil-mist-free, making them the best pal in laboratory.

### ► Quieter & more stable

Motor-driven transmission and vibration-proof assembly, making Rocker compressors much quieter and much stabler.

### ► Overheat protection

All Rocker compressors are equipped with overheat protection device that shuts down when overheated and resumes once cooled down, making the pumps safer and long-lasting.

### ► Air outlet filter

Every Rocker compressor is equipped with filter cartridge at the air outlet to filter particles and moisture to supply clean air.

### ► International certification

CE certification  
CSA certification  
IP 30 protection  
RoHS certification

### ► Two-Year Warranty

## Applications

- Aeration
- Pressure filtration

## Ordering Information

### 167320-11(22)

Rocker 320, Compressor  
AC110V, 60Hz (AC220V, 50Hz)

### 167420-11(22)

Rocker 420, Compressor  
AC110V, 60Hz (AC220V, 50Hz)



Rocker 320



Rocker 420

## Specification

Model	Rocker 320		Rocker 420	
PERFORMANCE DATA				
Max. pressure (w/o moisture trap)	90 psig		90 psig	
Max. pressure (w/ moisture trap)	80 psig		80 psig	
Max. flow rate	14 L/min @60Hz	12 L/min @50Hz	27 L/min @60Hz	25 L/min @50Hz
Hose barb	ID8		ID8	
Noise level	50 dB		52 dB	
Net weight	4.1 Kg		5.2 Kg	
Dimension (LxWxH)	27 x 14 x 17 cm		30 x 14 x 17 cm	
ELECTRICAL DATA				
Voltage	110V	220V	110V	220V
Frequency	50 / 60Hz		50 / 60Hz	
Max. power	90W (1/8 HP)	85W (1/8 HP)	147W (1/6 HP)	125W (1/6 HP)
Max. current	0.9A	0.4A	1.5A	0.6A
Safety	Overheat Protection		Overheat Protection	
OTHER DATA				
Equipped with	Pressure regulator Filter cartridge		Pressure regulator Filter cartridge	