

Pump Accessories



Hose Barb

Rocker Series Pump

Cat.No.:	167300-20	167300-21	167300-63	167300-23	167300-24	167300-25	167800-25
Barb :	3/8 inch	5/16 inch	1/4 inch	3/16 inch	1/8 inch	5/16 inch	3/8 inch
ID	ID10	ID8	ID6	ID5	ID3	ID8	ID10
Thread :	1/4 PT	1/8 PT	1/8 PT	1/8 PT	1/8 PT	1/8 PT	1/4 PT

Chemical Resistant Hose Barb

Chemker Series Pump

Cat.No.:	169300-21	169300-20
Barb :	5/16 inch	3/8 inch
ID	ID8	ID10
Thread :	1/4 NPS	1/4 NPS



Hose Barb

Tanker Series Pump

Cat.No.:	187100-32	187100-31	187100-34	167301-06	169411-06	189300-06
Barb :	1/2 inch	3/8 inch	1/4 inch	For Rocker	For Chemker	For Rocker
ID	ID12	ID10	ID6	300/400/410/430	300/400/410/411	300C/400C/410C
Thread :	1/4 SAE	1/4 SAE	1/4 SAE	167801-06	169611-06	189900-06
				For Rocker	For Chemker	For Rocker
				900/910	600/610/611	900C/910C



Pressure Regulator

167320-06	167440-06
For Rocker	For Rocker
320/420/430	440

Hoop Holder

167300-37	189300-37
For Rocker	For Rocker
300/400/410	300/400/410
300C/400C/410C	300C/400C/410C
Max. 118 mm	Max. 103 mm

Filter Cartridge

167301-07	167801-07
For Rocker	For Rocker
300/400/410/430	900/910
(with Prote™ overflow protection)	(without Prote™ overflow protection)



Check Valve

187100-33
For Tanker
130/150/215/230
ID10
1/4 SAE

Air Outlet (Copper Muffler)

167300-42
For Rocker
300/400/410/430
320/420

Silicone Tube

167600-42
For Rocker
900/901/910/911

Silicone Tube

167200-38	199100-68	180300-68
Ø 6 x 12 mm	Ø 4.2 x 8 mm	Ø 8 x 14 mm
Length: 100 cm	Length: 200 cm	Length: 200 cm

Universal Foot Switch

167200-41	167200-41-1
Momentary action	Maintained action

* PT (pipe thread) = BSPT