

121°C / 1bar / 20 min



# Filtration Apparatus Selection Guide

## SELECT THE BEST COMPONENTS FOR YOUR NEEDS


PPSU

### MF Pro, Magnetic Filter Holder

Model: MF 3 / 3a / 5 / 5a Pro

Feature: Magnetic seal

Page 35 ▶



Compatible with various bottles and flasks for direct filtrate collection.

PES

### MF, Magnetic Filter Holder

Model: MF 3 / 3a / 5 / 5a

Feature: Magnetic seal

Page 36 ▶



### LF, Spin-Lock Filter Holder

Model: LF 3 / 3a / 5 / 5a

Feature: Clamp-free lock

Page 37 ▶



### Filtration Apparatus

Model: LF 30 / MF 31

Feature: PES filtration set

Page 39 ▶



\* PPSU : Polyphenylsulfone, PES : Polyethersulfone



# Filtration Apparatus Selection Guide

SS 316

## SF, Spin-Lock Filter Holder

Model: SF 1 / 3 / 5  
Feature: Clamp-free lock



Page 38 ▶



## SFC, Continuous Filter Holder

Model: SF 142  
Feature: Continuous feed

Page 50 ▶



## Filtration Apparatus

Model: SF 10 / SF 30 / SF 11 / SF 52  
Feature: SS316 filtration set

Page 40 ▶



Page 41 ▶



Borosilicate Glass

## VFC, Continuous Filter Holder

Model: VFC 47  
Feature: Continuous feed

Page 49 ▶



## Continuous Filter Holder

Model: Stream 3 / 10  
Feature: Continuous feed

Page 48 ▶



## 25 mm Filtration Apparatus

Model: VF 5 / 8  
Feature: Small-volume sample filtration

Page 42 ▶



## 47 mm Filtration Apparatus

Model: VF 3  
Feature: All-glass construction

Page 45 ▶



## 47 mm Filtration Apparatus

Model: VF 6 / 7  
Feature: General filtration with secure fixing

Page 43 ▶



## 47 mm Filtration Apparatus

Model: VF 12 / 14  
Feature: Filter and store in one step

Page 44 ▶



## 47 mm Filtration Apparatus

Model: VF 20 / 21  
Feature: Filter in evaporating flask

Page 46 ▶



## 90 mm Filtration Apparatus

Model: VF 15 / VF 16 / VF 17  
Feature: Filter and store in one step

Page 47 ▶

