

# SF, Stainless Steel Filter Holder



## Features

- ▶ Clamp-free spin-lock design.
- ▶ Precise laser-marking ensures  $\leq 2.5\%$  volume error at max. graduation.
- ▶ Membrane guide at the support base helps position membrane easily.
- ▶ Forceps access points on support base help retrieve membrane without damage.
- ▶ SS316 funnel is shatterproof, durable and can be fast sterilized by flame.
- ▶ Two-Year Warranty

*Recommended for microbiological applications that require fast flame sterilization.*

## Specification

Material	SF 1	SF 3	SF 5
Funnel	SS316	SS316	SS316
Support base	SS316	SS316	SS316
Membrane support	SS316	SS316	SS316
Funnel lid	SS316	SS316	SS316

### Specification

Funnel capacity	100 mL	300 mL	500 mL
Filter diameter	47 mm	47 mm	47 mm
Effective filtration area	11.6 cm <sup>2</sup>	11.6 cm <sup>2</sup>	11.6 cm <sup>2</sup>
Graduated increment	50 / 100 mL	100 / 250 / 300 mL	250 / 500 mL

## Ordering Information

SF 1, Stainless Steel Filter Holder (100 mL) 180100-10

1	Stainless steel funnel, 100 mL	180100-01
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16

SF 3, Stainless Steel Filter Holder (300 mL) 180100-30

1	Stainless steel funnel, 300 mL	180100-03
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16

SF 5, Stainless Steel Filter Holder (500 mL) 180100-50

1	Stainless steel funnel, 500 mL	180100-05
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16

### Optional

Stainless steel funnel lid, 100 mL	180100-11
Stainless steel funnel lid, 300 mL / 500 mL	180100-13

