

Stainless Steel Filtration Apparatus

121°C / 1bar / 20 min
Autoclavable

2
YEAR
WARRANTY



SF 10



SF 30

Features

- ▶ Handy, clamp-free SF filter holder with spin-lock design.
- ▶ SS316 funnel is shatterproof, durable and can be fast sterilized by flame.
- ▶ Built-in overflow protection makes filtration easier and care-free.
- ▶ Two-Year Warranty (SF Filter Holder)

Applications


General filtration, such as medium purification, microbiological test etc.

Specification

Material	SF 10	SF 30
Funnel	SS316	SS316
Support base	SS316	SS316
Membrane support	SS316	SS316
Receiver bottle	PC	PC

Specification

Funnel capacity	100 mL	300 mL
Bottle capacity	1200 mL	1200 mL
Filter diameter	47 mm	47 mm
Effective filtration area	11.6 cm ²	11.6 cm ²
Hose barb	ID8	ID8
Overflow protection	Yes	Yes



SF 10

SF 30

Ordering Information

SF 10, Stainless Steel Filtration Set 197010-32

1	Stainless steel funnel, 100 mL	180100-01
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.16)	167110-17
5	PC cover set	197000-32-PC
6	PC bottle, 1200 mL	197000-31-PC

SF 30, Stainless Steel Filtration Set 197010-33

1	Stainless steel funnel, 300 mL	180100-03
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.16)	167110-17
5	PC cover set	197000-32-PC
6	PC bottle, 1200 mL	197000-31-PC