

Maximize Your MultiVac Versatility!

A

Step 1. Choose the manifold



MultiVac 300 - MB
3 Places Stainless Steel Manifold
(180300-02)



MultiVac 600 - MB
6 Places Stainless Steel Manifold
(180600-02)



Step 2. Choose the funnel



ROCKER, 100 mL
SF stainless steel funnel
(180100-01)



Pall, 100 mL / 250 mL
Sentino® disposable plastic funnel
(4870 / 4871)



ROCKER, 300 mL
SF stainless steel funnel
(180100-03)



Stainless steel funnel lid
for SF funnel, 100 mL
(180100-11)



ROCKER, 500 mL
SF stainless steel funnel
(180100-05)



Stainless steel funnel lid
for SF funnel, 300 mL / 500 mL
(180100-13)

B

Step 1. Choose the manifold



MultiVac 300 - MS
3 Places Stainless Steel Manifold
(180300-01)



MultiVac 600 - MS
6 Places Stainless Steel Manifold
(180600-01)



Step 2. Choose the filter holder



ROCKER, 300 mL
MF 3 / MF 3 Pro,
Magnetic Filter Holder
(200300-01 / 200300-01-P)



ROCKER, 500 mL
MF 5 / MF 5 Pro,
Magnetic Filter Holder
(200500-01 / 200500-01-P)



ROCKER, 300 mL
Glass Filter Holder
(167120-31)

* MB adaptors are compatible with Pall, Sentino® disposable funnels (4870 / 4871).

C

Step 1. Choose the manifold




MultiVac 300
3 Places Stainless Steel Manifold Base
(180300-00)



MultiVac 600
6 Places Stainless Steel Manifold Base
(180600-00)

Customized Design

We can offer customized adaptors for your funnel



Step 2. Choose the disposable funnel and adaptor



Sartorius, 100 mL / 250 mL
Microsart® disposable plastic funnel
(16A07-10-N / 16A07-25-N)



Millipore, 100 mL / 250 mL
Microfil® disposable plastic funnel
(MIHAWG100 / MIHAWG250)



Silicone tube
Ø 4.2 x 8 mm, 200 cm
(199100-68)

+



SS316 Adaptor for Sartorius
disposable plastic funnel
(190100-74)

+



SS316 Adaptor for Millipore
disposable plastic funnel
(190100-73)

+



SS316 Hose barb
for silicone tube (ID4 ~ ID8)
(190100-76)



* SS316 adaptor is only available for above-mentioned disposable funnels.