


# Suction Accessories



## PP Heavy-Duty Vacuum Bottle

| Specification | PP vacuum bottle, 4000 mL                                                         | PP vacuum bottle with drain, 4000 mL                                              | PP vacuum bottle with quick-coupling connector, 4000 mL                            | PP vacuum bottle, 4000 mL                                                           |
|---------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|               |  |  |  |  |

|                     |           |                                |                          |           |
|---------------------|-----------|--------------------------------|--------------------------|-----------|
| Catalog number      | 197200-47 | 197200-48                      | 197200-53                | 197200-46 |
| Connector           | Hose barb | Hose barb                      | Quick-coupling connector | -         |
| Graduated increment | 500 mL    | 500 mL                         | 500 mL                   | 500 mL    |
| Overflow protection | Yes       | Yes                            | Yes                      | -         |
| Hose barb           | ID4 ~ ID8 | ID4 ~ ID8                      | ID6                      | -         |
| Drain device        | -         | Yes (Rubber, NOT autoclavable) | -                        | -         |
| Autoclavable        | Yes       | Yes                            | Yes                      | Yes       |



Refer to P.79  
for full PP bottle series.

\* PP : Polypropylene



# Suction Accessories



## PP Vacuum Bottle

| Specification       | PP vacuum bottle, 1000 mL | PP vacuum bottle, 2000 mL |
|---------------------|---------------------------|---------------------------|
| Catalog number      | 197100-10                 | 197202-47                 |
| Graduated increment | 250 mL                    | 500 mL                    |
| Overflow protection | -                         | Yes                       |
| Hose barb           | ID6                       | ID4 ~ ID8                 |
| Autoclavable        | Yes                       | Yes                       |



## PES Vacuum Bottle

| Specification            | PES vacuum bottle, 1200 mL | PES vacuum bottle, 1200 mL |
|--------------------------|----------------------------|----------------------------|
| Catalog number           | 197000-11                  | 197000-12                  |
| Graduated increment      | 100 mL                     | 100 mL                     |
| Overflow protection      | Yes                        | Yes                        |
| Cap with suction fitting | Yes                        | -                          |
| Hose barb                | ID8                        | ID4 ~ ID8                  |
| Autoclavable             | Yes                        | Yes                        |



## Quick Connect Fitting

| Specification  | 1/4 inch Male fitting | 5/16 inch Female fitting |
|----------------|-----------------------|--------------------------|
| Catalog number | 197201-38             | 197201-39                |
| Hose barb      | ID6                   | ID8                      |
| Material       | POM, FKM              | POM, SS316, FKM          |
| Autoclavable   | Yes                   | Yes                      |



## Laboratory Glass Bottle with GL45 Suction Cap

| Specification       | Glass bottle, 500 mL | Glass bottle, 1000 mL | Glass bottle, 3500 mL |
|---------------------|----------------------|-----------------------|-----------------------|
| Catalog number      | 197000-60-500        | 197000-60-1000        | 197000-60-3500        |
| Cap                 | GL45 suction cap     | GL45 suction cap      | GL45 suction cap      |
| Graduated increment | 100 mL               | 100 mL                | 500 mL                |
| Hose barb           | ID6                  | ID6                   | ID6                   |
| Autoclavable        | Yes                  | Yes                   | Yes                   |

\* PES : Polyethersulfone, PP : Polypropylene