


Suction Accessories

PP Heavy-Duty Vacuum Bottle

Specification	4000 ml, PP vacuum bottle	4000 ml, PP vacuum bottle with drain	4000 ml, PP vacuum bottle with quick-coupling connector
---------------	---------------------------	--------------------------------------	---



Catalog number	197200-47	197200-48	197200-53
Connector	Hose barb	Hose barb	Quick-coupling connector
Graduated increment	500 ml	500 ml	500 ml
Overow protection	Yes	Yes	Yes
Hose barb	ID4 ~ ID8	ID4 ~ ID8	ID6
Drain device	-	Yes	-
Autoclavable	Yes	Yes	Yes



* PP : Polypropylene