

MF, Magnetic Filter Holder

One to Many

Easy

Reusable



Features



One to Many

- 1 ♦ Compatible with various kinds of bottle, vacuum flask with or without side arm to collect filtrate directly.



One-handed & Easy Operation

- 2 ♦ Strong magnetic seal allows one-handed operation.
♦ Membrane guide helps position membrane easily.
♦ Forceps access points help retrieve membrane easily.



Sturdy & Reusable

- 3 ♦ PPSU (Polyphenylsulfone) and PES (Polyethersulfone) construction provides good chemical resistance.
♦ High temperature resistance withstands autoclave cycles.

Applications

- ♦ General filtration
- ♦ Water testing



Certification

RoHS

Compatible with various kinds of bottle and vacuum flask



Specification

Model	MF 3 / 3a (Pro)	MF 5 / 5a (Pro)
Funnel capacity	300 ml	500 ml
Effective filtration area	13.1 cm ²	13.1 cm ²
Filter diameter	47 mm	47 mm
Graduated increment	50 ml	50 ml
Hose barb	ID4 ~ ID8	ID4 ~ ID8

Material	MF 3/ 3a / 5 / 5a Pro	MF 3/ 3a / 5 / 5a
Funnel	PPSU	PES
Support base	PPSU	PES
Membrane support	PPSU	PES
Hose barb	PPSU	PES
Funnel lid	PP	PP
Cap of lid	Silicone	Silicone
Syringe filter	PTFE	PTFE



MF 3 Pro, 300ml Magnetic Filter Holder with lid kit



MF 5 Pro, 500ml Magnetic Filter Holder with lid kit



MF 3a Pro, 300ml Magnetic Filter Holder (short stem) with lid kit



MF 5a Pro, 500ml Magnetic Filter Holder (short stem) with lid kit



200300-01-P	MF 3 Pro, Magnetic Filter Holder (300 ml)
200300-01	MF 3, Magnetic Filter Holder (300 ml)
200500-01-P	MF 5 Pro, Magnetic Filter Holder (500 ml)
200500-01	MF 5, Magnetic Filter Holder (500 ml)
200310-01-P	MF 3a Pro, Magnetic Filter Holder (300 ml)
200310-01	MF 3a, Magnetic Filter Holder (300 ml)
200510-01-P	MF 5a Pro, Magnetic Filter Holder (500 ml)
200510-01	MF 5a, Magnetic Filter Holder (500 ml)
200300-36	PP Lid Kit
167110-16	Silicone Stopper (No. 8)
167100-52	PALL, GN-6 Membrane, grid, sterilized, 0.45 um
167100-20	PALL, A/E Glass fiber Filter (1 um), 100/PK

Rocker Scientific Co., Ltd.

Headquarter

19F.-6, No. 206, Guanghua 1st Rd., Lingya Dist.,
Kaohsiung City 802040, Taiwan (R.O.C.)
T/ +886-7-227-3358 F/ +886-7-227-3716

Taipei office

11F., No. 402, Sec. 1, Ren'ai Rd., Linkou Dist.,
New Taipei City 244014, Taiwan (R.O.C.)
T/ +886-2-2603-3311 F/ +886-2-2603-6622

Get more information for Rocker

www.rocker.com.tw
export@rocker.com.tw

