

Oil Free Vacuum Pump / Compressor



Features

- ▶ **2-in-1 machine**
Rocker 430 works as both vacuum pump and compressor. Equipped with gauges and regulators for easy monitoring and adjusting the vacuum level or delivery pressure.
- ▶ **Oil-free & care-free**
Rocker pumps are driven by piston that needs no lubricant, making them maintenance-free and oil-mist-free.
- ▶ **Prote™ overflow protection**
Equipped with Prote™ overflow protection to prevent liquid from being sucked into the pump.
- ▶ **Overheat protection**
All Rocker pumps are equipped with overheat protection device that shuts down when overheated and resumes once cooled down, making the pumps safer and long-lasting.
- ▶ **Moisture trap with filter cartridge**
Filter cartridge filters particles and moisture to ensure air cleanness of both inlet and supply.
- ▶ **International certification**
CE certification
IP 30 protection
RoHS certification
- ▶ **Two-Year Warranty**

Applications

- ▶ Vacuum filtration
- ▶ Microbiological test
- ▶ Suspended solid test
- ▶ Air sampling
- ▶ Pressure filtration

Ordering Information

167430-11(22)

Rocker 430, Vacuum Pump
AC110V, 60Hz (AC220V, 50Hz)



Specification

Model	Rocker 430	
PERFORMANCE DATA		
Max. vacuum	145 mbar abs.	
Max. pressure	80 psig	
Max. flow rate	34 L/min @60Hz	28 L/min @50Hz
Hose barb	ID8	
Noise level	52 dB	
Net weight	5.5 Kg	
Dimension (LxWxH)	25 x 24 x 20 cm	
ELECTRICAL DATA		
Voltage	110V	220V
Frequency	50 / 60Hz	
Max. power	147W (1/6 HP)	125W (1/6 HP)
Max. current	1.5A	0.6A
Safety	Overheat Protection	
OTHER DATA		
Equipped with	Vacuum / Pressure regulator Overflow protection / Filter cartridge	