

# Vacuum Filtration System



## Features

### ► For large number of samples

MultiVac 310 / 610 - MS - T is an economical kit for filtration with large number samples.

### ► Light and handy

This aluminum-made manifold is light, cost-effective and autoclavable. Each component can be easily removed for cleaning and autoclaving.

### ► Magnetic filter holder is included

Magnetic filter holder allows for one-handed operation, and can prevent the membrane from being twisted or torn.

### ► 4000 mL PP vacuum bottle with drain

The large capacity of 4000 mL vacuum bottle with drain offers convenience of collecting and draining liquid waste.

### ► International certification

CE certification

CSA certification (Rocker 400)

IP 30 protection (Rocker 400)

RoHS certification (Rocker 400)

### ► Two-Year Warranty (Rocker 400)

### ► Two-Year Warranty (MultiVac 310 / 610 - MS)

## MultiVac 310 - MS - T



## MultiVac 610 - MS - T



## Applications

- Suspended solid test
- Microbiological test
- Vacuum filtration

## Ordering Information

**180311-11(22)-T**

**180611-11(22)-T**

MultiVac 310 / 610 - MS - T, Vacuum Filtration System, complete set includes :

\* Rocker 400, Oil Free Vacuum Pump

\* MultiVac 310 / 610 - MS, 3/6 Places

Aluminum Manifold

\* MF 3, Magnetic Filter Holder, 300 mL (3/6 pcs)

\* PP Vacuum Bottle with drain, 4000 mL

\* Silicone Tube (2 pcs)

## Specification

Model	MultiVac 310 - MS - T		MultiVac 610 - MS - T	
Voltage	110V	220V	110V	220V
Max. power	100W	80W	100W	80W
Max. vacuum	99 mbar abs.			
Max. flow rate	37 L/min @60Hz	34 L/min @50Hz	37 L/min @60Hz	34 L/min @50Hz
Vacuum regulator	Yes			
Overflow protection	Yes			
Noise level	52 dB			
OTHER DATA				
Number of branch	3		6	
Filter holder	MF 3, 300 mL x 3		MF 3, 300 mL x 6	
Filter diameter	47 mm		47 mm	
Material contact with filtrate	Aluminum		Aluminum	