

Lab Burner



Features


- ▶ **Reliable & professional**
 - (1) Metal-made base ensures steadiness and spill-proof.
 - (2) Compliant to CE safety regulation for gas-fired appliance.
- ▶ **Stable flame**
 - (1) Superior gas inlet design ensures sufficient butane feed and stable flame outputs.
 - (2) Burns up to 1300°C for max. 60 minutes.
- ▶ **Easy yet safe**
 - (1) Triggered by Piezo Igniter, providing instant flame at one press.
 - (2) Additional safety lock provides additional security. (Dragon 100)
- ▶ **One-Year Warranty**

Applications

- ▶ Sterilizing
- ▶ Welding
- ▶ Heating
- ▶ Cooking

Ordering Information

- 177100-00**
Dragon 100, Torch Burner
- 177200-00**
Dragon 200, Lab Burner
- 177220-00**
Dragon 220, Lab Burner
- 177200-T**
Dragon 200 - T, Lab Burner Set includes:
 - * Dragon 200, Lab Burner
 - * Lab Burner Stand
 - * Ceramic Net Mesh (14 x 14 cm)
- 177200-10**
Lab Burner Stand
- 177200-12**
Ceramic Net Mesh (14 x 14 cm)



Dragon 100



Dragon 200 - T



Dragon 200



Dragon 220

Specification

Model	Dragon 100	Dragon 200	Dragon 200 - T	Dragon 220
PERFORMANCE DATA				
Max. temperature	1300 °C	1300 °C	1300 °C	1300 °C
Fuel tank capacity	40 g	35 g	35 g	35 g
Working time	60 min	60 min	60 min	60 min
Safe switch	Yes	-	-	-
Net weight	380 g	410 g	1235 g	425 g
Dimension (LxWxH)	8 x 14 x 16 cm	8 x 12 x 9 cm	15 x 17 x 13 cm	8 x 12 x 9 cm
OTHER DATA				
Equipped with	-	-	Lab burner stand Ceramic net mesh	-