

Electric Sterilizer

Features

► **Fuel-free**

Dragon 300 / 320 are electric burners for rapid sterilization. By generating heat from electricity, Dragon 300 / 320 is fuel-free and free of the risks of cross-contamination that may be caused by open flame.

► **Fast & stable heating**

Max. temperature is up to 850°C. Powerful heating elements ensure fast and stable heating.

► **Wide entrance**

With much wider entrance, Dragon 320 can hold objects with larger head diameter up to 37 mm.

► **User-friendly & easy**

Weighted base, loop rack and adjustable chamber angle (from 45° to 75°) make Dragon 300 / 320 handy and easy to use inside laminar flow hoods or anaerobic chambers.

► **International certification**

CE certification

► **One-Year Warranty**

Applications

- Culture tube sterilization
- Inoculating loops and needles sterilization

Ordering Information

177300-11(22)

Dragon 300, Electric Sterilizer
AC110V, 60Hz (AC220V, 50Hz)

177320-11(22)

Dragon 320, Electric Sterilizer
AC110V, 60Hz (AC220V, 50Hz)

Dragon 300



Dragon 320



Specification

| Model | Dragon 300 | | Dragon 320 | |
|-------------------------|----------------|------|----------------|------|
| PERFORMANCE DATA | | | | |
| Max. temperature | 850 °C | | | |
| Standby temperature | 600 °C | | | |
| Chamber diameter | 16 mm | | 37 mm | |
| Net weight | 1.1 Kg | | 1.1 Kg | |
| Dimension (LxWxH) | 17 x 9 x 21 cm | | 19 x 9 x 22 cm | |
| ELECTRICAL DATA | | | | |
| Voltage | 110V | 220V | 110V | 220V |
| Frequency | 50 / 60Hz | | 50 / 60Hz | |
| Max. power | 180W | | 180W | |
| Safety | Fuse | | Fuse | |