

Air Supply System



Features

► **Compact & care-free**

Rocker oil-free compressors are compact in size, light in weight and oil-mist-free, making them the best pal in laboratory.

► **Quieter & more stable**

Motor-driven transmission and vibration-proof assembly, making Rocker compressors much quieter and stabler.

► **Overheat protection**

All Rocker compressors are equipped with overheat protection device that shuts down when overheated and resumes once cooled down, making the pumps safer and long-lasting.

► **Air outlet filter**

Every Rocker series compressor is equipped with filter cartridge at the air outlet to filter particles and moisture to supply clean air.

► **International certification**

CE certification
CSA certification
IP 30 protection
RoHS certification

► **Two-Year Warranty**

Applications

- Fermenter
- Bioreactor
- Gas chromatograph

Ordering Information

167440-11(22)

Rocker 440, Air Supply System
AC110V, 60Hz (AC220V, 50Hz)

167304-31

Auto Stop Switch, 60~80 psig



Rocker 440

*The capacity of SS304 tank is 2.5 L

Specification

Model		Rocker 440	
PERFORMANCE DATA			
Pressure range	40~60 psig		
Max. flow rate (w/o tank)	27 L/min @60Hz	25 L/min @50Hz	
Max. flow rate (w/ tank)	20 psig	40 L/min @60Hz	40 L/min @50Hz
	40 psig	52 L/min @60Hz	50 L/min @50Hz
	60 psig	65 L/min @60Hz	60 L/min @50Hz
Hose barb	ID8		
Noise level	52 dB		
Net weight	8.4 Kg		
Dimension (LxWxH)	29 x 31 x 25 cm		
ELECTRICAL DATA			
Voltage	110V	220V	
Frequency	50 / 60Hz		
Max. power	147W (1/6 HP)	125W (1/6 HP)	
Max. current	1.5A	0.6A	
Safety	Overheat Protection		
OTHER DATA			
Equipped with	Pressure regulator Filter cartridge		

*The pressure range can be regulated by replacing the auto stop switch.