

SF, Stainless Steel Filter Holder



Features

- ▶ Clamp-free spin-lock design.
- ▶ Precise laser-marking ensures high volume accuracy.
- ▶ Membrane guide at the support base helps position membrane easily.
- ▶ Forceps access points on support base help retrieve membrane without damage.
- ▶ SS316 funnel is shatterproof, durable and can be fast sterilized by flame.
- ▶ Two-Year Warranty

Recommended for microbiological applications that required fast flame sterilization.

Specifications

Material	SF 1	SF 3	SF 5
Funnel	SS316	SS316	SS316
Support base	SS316	SS316	SS316
Membrane support	SS316	SS316	SS316
Funnel lid	SS316	SS316	SS316

Specifications

Funnel capacity	100 mL	300 mL	500 mL
Filter diameter	47 mm	47 mm	47 mm
Effective filtration area	11.6 cm ²	11.6 cm ²	11.6 cm ²
Graduated increment	50 / 100 mL	100 / 250 / 300 mL	250 / 500 mL

Ordering Information

SF 1, Stainless Steel Filter Holder (100 mL) 180100-10

1	Stainless steel funnel, 100 mL	180100-01
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16

SF 3, Stainless Steel Filter Holder (300 mL) 180100-30

1	Stainless steel funnel, 300 mL	180100-03
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16

SF 5, Stainless Steel Filter Holder (500 mL) 180100-50

1	Stainless steel funnel, 500 mL	180100-05
2	Stainless steel support screen	180100-31
3	Stainless steel support base	180100-00
4	Silicone stopper (No.8)	167110-16

SF, Stainless Steel Funnel Lid

Stainless steel funnel lid, 100 mL	180100-11
Stainless steel funnel lid, 300 mL / 500 mL	180100-13

