

Suction Accessories




PP Heavy-Duty Vacuum Bottle

Specification	PP vacuum bottle, 4000 mL	PP vacuum bottle with drain, 4000 mL	PP vacuum bottle with quick-coupling connector, 4000 mL
---------------	---------------------------	--------------------------------------	---



Catalog number	197200-47	197200-48	197200-53
Connector	Hose barb	Hose barb	Quick-coupling connector
Graduated increment	500 mL	500 mL	500 mL
Overow protection	Yes	Yes	Yes
Hose barb	ID4 ~ ID8	ID4 ~ ID8	ID6
Drain device	-	Yes	-
Autoclavable	Yes	Yes	Yes



* PP : Polypropylene