

Suction Accessories



PES Vacuum Bottle

Specification	PES vacuum bottle, 1200 mL	PES vacuum bottle, 1200 mL
Catalog number	197000-11	197000-12
Graduated increment	100 mL	100 mL
Overflow protection	Yes	Yes
Cap with suction fitting	Yes	-
Hose barb	ID8	ID4 ~ ID8
Autoclavable	Yes	Yes



PP Vacuum Bottle

Specification	PP vacuum bottle, 1000 mL	PP vacuum bottle, 2000 mL
Catalog number	197100-10	196100-0406
Graduated increment	250 mL	-
Overflow protection	-	-
Hose barb	ID6	ID6
Autoclavable	Yes	Yes



Quick Connect Fitting

Specification	1/4 inch Male Fitting	1/4 inch Female Fitting
Catalog number	197201-38	197201-39
Hose barb	ID6	ID6
Material	POM, FKM	POM, SS316, FKM
Autoclavable	Yes	Yes



GL45 Laboratory Glass Bottle

Specification	Glass bottle, 500 mL	Glass bottle, 1000 mL	Glass bottle, 3500 mL
Catalog number	197000-60-500	197000-60-1000	197000-60-3500
Cap	GL45 suction cap	GL45 suction cap	GL45 suction cap
Graduated increment	100 mL	100 mL	500 mL
Hose barb	ID6	ID6	ID6
Autoclavable	Yes	Yes	Yes

* PES : Polyethersulfone, PP : Polypropylene