

Filtration Accessories



Membrane Filters

Specification	Glass fiber filter	MCE membrane
Catalog number	167100-20	167100-52
Brand	Pall	Pall
Filter diameter	47 mm	47 mm
Pore size	1 µm	0.45 µm
Material of filter	Glass fiber	MCE
Package	100/PK	200/PK



GL45 Filtration Adaptor

Specification	GL45 filtration adaptor
Catalog number	197000-65
Material of cap	PP
Material of gasket	PTFE
Stopper	No.8 stopper
Hose barb	ID4 ~ ID8
Autoclavable	Yes



GL45 Suction Cap NEW

Specification	GL45 suction cap	SC45 Suction Cap	SC45 Suction Cap Kit
Catalog number	197000-60	197000-61-P2N	197000-61-P2
Material of cap	PP	PP	PP
Material of gasket / O-ring	Silicone	EPDM	EPDM
Overflow protection	-	-	Yes
Hose barb	ID6	ID4 ~ ID8	ID4 ~ ID8
Autoclavable	Yes	Yes	Yes

Application for GL45 Suction Cap

On waste bottle

Collecting liquid waste by being connected to vacuum pumps and filtration accessories.



On moisture trap bottle

Attaching to GL45 bottle filled with desiccant to prevent moisture from entering the vacuum pump.



On safety bottle

Working as a safety bottle to prevent liquid from being sucked into pump.



* Please refer to P.82 to select GL45 Laboratory Glass Bottle