

# Electric Sterilizer

## Features

### ► Fuel-free

Dragon 300 / 320 are electric burners for rapid sterilization. By generating heat from electricity, Dragon 300 / 320 is fuel-free and free of the risks of cross-contamination that may caused by open flame.

### ► Fast & stable heating

Max. temperature is up to 850°C. Powerful heating elements ensure fast and stable heating.

### ► Wide entrance

With much wider entrance, Dragon 320 can hold objects with larger head diameter up to 37 mm.

### ► User-friendly & easy

Weighted base, loop rack and adjustable chamber angle (from 45° to 75°) make Dragon 300 / 320 handy and easy to use inside laminar flow hoods or anaerobic chambers.

### ► International certification

CE certification

### ► One-Year Warranty

## Applications

- Culture tube sterilization
- Inoculating loops and needles sterilization

## Ordering Information

### 177300-11(22)

Dragon 300, Electric Sterilizer  
AC110V, 60Hz (AC220V, 50Hz)

### 177320-11(22)

Dragon 320, Electric Sterilizer  
AC110V, 60Hz (AC220V, 50Hz)

### Dragon 300



### Dragon 320



## Specification

Model	Dragon 300		Dragon 320	
PERFORMANCE DATA				
Max. temperature	850 °C			
Standby temperature	600 °C			
Chamber diameter	16 mm		37 mm	
Net weight	1.1 Kg		1.1 Kg	
Dimension (LxWxH)	17 x 9 x 21 cm		19 x 9 x 22 cm	
ELECTRICAL DATA				
Voltage	110V	220V	110V	220V
Frequency	50 / 60Hz		50 / 60Hz	
Max. power	180W		180W	
Safety	Fuse		Fuse	