



PP Heavy-Duty Vacuum Bottle



PP Heavy-Duty Vacuum Bottle, 4000 mL

Specification				
Catalog number	197200-47	197200-48	197200-53	197200-46
Cap size	B83	B83	B83	B83
Connector	Hose barb	Hose barb	Quick-coupling connector	-
Graduated increment	500 mL	500 mL	500 mL	500 mL
Overflow protection	Yes	Yes	Yes	-
Hose barb	ID4 ~ ID8	ID4 ~ ID8	ID6	-
Drain device	-	Yes	-	-



PP Heavy-Duty Vacuum Bottle, 2000 mL

Specification			
Catalog number	197202-47	197202-53	197202-46
Cap size	B83	B83	B83
Connector	Hose barb	Quick-coupling connector	-
Graduated increment	500 mL	500 mL	500 mL
Overflow protection	Yes	Yes	-
Hose barb	ID4 ~ ID8	ID6	-



PP Heavy-Duty Vacuum Bottle, 1000 mL

Specification			
Catalog number	197100-10	197100-10SC	197100-10SCN
Cap size	GL45	GL45	GL45
Connector	Hose barb	Hose barb	Hose barb
Graduated increment	250 mL	250 mL	250 mL
Overflow protection	-	Yes	-
Hose barb	ID6	ID4 ~ ID8	ID4 ~ ID8