

Suction Kit



121°C / 1bar / 20 min

Autoclavable

Features

- ▶ **Quick connection**
Easy installation between handle and adaptors.
- ▶ **One-finger tip ejection**
Ergonomic ejector makes tip removal easy and effortless.
- ▶ **Easy lock switch**
Keep continuous suction with a single click to save efforts during repetitive procedures.
- ▶ **Multi-funtional storage rack**
Accommodating all accessories in this rack to keep a tidy working environment. Built-in liquid tray collects dripping liquid from the adaptors and can be cleaned easily.
- ▶ **Autoclavable**
The complete kit except the storage rack are made of autoclavable materials including POM and SS304.

Applications

- ▶ Cell culture aspiration
- ▶ Laboratory waste suction

Ordering Information

199100-00

BioDolphin Suction Kit (Light Blue)

199200-00

Plus Suction Kit (Light Blue)

199110-00

BioDolphin Suction Kit (Deep Blue)

199210-00


Plus Suction Kit (Deep Blue)

199130-00

BioDolphin Suction Kit (Green)

199230-00

Plus Suction Kit (Green)



Standard accessory :

- | | |
|---|---|
| 1. 199100-30 Handle | 6. 199100-49 8-channel SS adaptor (40 mm) |
| 2. 199100-46 1-channel tip adaptor (50 mm) | 7. 199100-75 1-channel tip adaptor with ejector |
| 3. 199100-76 1-channel tip adaptor (150 mm) | 8. 199100-50 8-channel tip adaptor with ejector |
| 4. 199100-47 1-channel SS adaptor (40 mm) | 9. 199100-61 Rack |
| 5. 199100-48 1-channel SS adaptor (80 mm) | |

Specification

Model	BioDolphin	Plus
Operation	Press to suck	
Hose barb	ID6	
Suction kit material	POM, SS304	
Rack material	ABS	-
Net weight	600 g	100 g
Dimension (LxWxH)	15 x 15 x 20 cm	15 x 8 x 6 cm
Suction kit	1. Handle 2. 1-channel tip adaptor (50 mm) 3. 1-channel tip adaptor (150 mm) 4. 1-channel SS adaptor (40 mm) 5. 1-channel SS adaptor (80 mm) 6. 8-channel SS adaptor (40 mm) 7. 1-channel tip adaptor with ejector 8. 8-channel tip adaptor with ejector 9. Rack	1. Handle 3. 1-channel tip adaptor (150 mm) 7. 1-channel tip adaptor with ejector 8. 8-channel tip adaptor with ejector

*The storage rack is NOT autoclavable.