



Vacuum Filtration System new evolutions of **HOOP**

Safe · Space-saving · Simple

Features



Safe

- 1 ♦ Hoop design prevents the negligence from tipping over the waste bottle.



Space-saving

- 2 ♦ Solve tubing stress, 50% space-saving than traditional DIY system.



Simple

- 3 ♦ No air pollution, maintenance free.
♦ Fix the waste bottle in 3 easy steps.

Applications

- ♦ Vacuum filtration
- ♦ Suspended solid test
- ♦ Solvent purification (Rocker 300C - VF 12)



Certification



Lift the Hoop



Place the filtration apparatus into Hoop



Press the Hoop down

Online video



Specification

Model	Rocker 300C - VF 12		Rocker 300 - SF 10 / SF 30 / LF 30		
Voltage / Freq.	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz	
Max.vacuum	99 mbar		99 mbar		
Max. flow rate	23 L/min	20 L/min	23 L/min	20 L/min	
Vacuum regulator	Yes		Yes		
Overflow protection / Filter cartridge	-		Yes		
Hose barb	ID8		ID8		
Model	VF 12		SF 10	SF 30	LF 30
Funnel capacity	300 ml		100 ml	300 ml	
Bottle capacity	1000 ml		1200 ml		
Filter diameter	47 mm		47 mm		



Rocker 300C - VF 12



Rocker 300 - SF 10



Rocker 300 - SF 30



Rocker 300 - LF 30



189300-37	Hoop holder for Rocker 300C - VF 12
167300-37	Hoop holder for Rocker 300 - LF 30 / LF 30 - SS / SF 30 / SF 10
189302-11(22)	Rocker 300C - VF12, Solvent purification system, AC110V, 60Hz (AC220V, 50Hz)
167305-11(22)	Rocker 300 - LF 30, Vacuum filtration system, AC110V, 60Hz (AC220V, 50Hz)
167309-11(22)	Rocker 300 - SF 30, Vacuum filtration system, AC110V, 60Hz (AC220V, 50Hz)
167306-11(22)	Rocker 300 - SF 10, Vacuum filtration system, AC110V, 60Hz (AC220V, 50Hz)
167200-12	VF 12, Glass filtration set
197010-30	LF 30, PES filtration set
197010-32/33	SF 10 / 30, Stainless steel filtration set
167300-11(22)	Rocker 300, Vacuum pump, AC110V, 60Hz (AC220V, 50Hz)
189300-11(22)	Rocker 300C, PTFE coated chemical resistant vacuum pump, AC110V, 60Hz (AC220V, 50Hz)
167200-38	Silicone tube, Ø 6 x 12 mm, 1 m
167100-20	PALL, Type A/ E Glass Fiber Discs, 1 µm, 100/PK
167100-52	PALL, GN-6 MCE membrane, grid, 0.45 µm, 200/pk

Rocker Scientific Co., Ltd.

Headquarter

 19F.-6, No.206, Guanghua 1st Rd., Lingya Dist.,
 Kaohsiung City 802, Taiwan (R.O.C.)
 T/ +886-7-227-3358 F/ +886-7-227-3716

Taipei office

 11F., No. 402, Sec. 1, Ren'ai Rd., Linkou Dist.,
 New Taipei City 244014, Taiwan (R.O.C.)
 T/ +886-2-2603-3311 F/ +886-2-2603-6622

Get more information for Rocker

www.rocker.com.tw
export@rocker.com.tw
